



TruenatTM

MTB

Chip-based Real Time PCR test for *Mycobacterium tuberculosis*

1. INTENDED USE


TruenatTM MTB (REF 601030005 / 601030020) is a chip-based Real Time Polymerase Chain Reaction (PCR) test for the quantitative detection and diagnosis of *Mycobacterium tuberculosis* (MTB) in human pulmonary (sputum/non-sputum) and EPTB specimen and aids in the diagnosis of infection with MTB. **TruenatTM MTB** runs on the **TruelabTM Uno** and **TruelabTM Uno Dx** Real Time micro PCR Analyzer.

2. INTRODUCTION

Tuberculosis (TB) is an infectious disease caused predominantly by the bacillus *Mycobacterium tuberculosis*. It typically affects the lungs (pulmonary TB) but can affect other sites as well (extra pulmonary TB). Tuberculosis (TB) is the second largest killer worldwide, after HIV and is the leading cause of death in HIV patients. Pulmonary TB spreads through air and is highly contagious. Over 80% of TB infections are pulmonary and if left untreated, a pulmonary TB patient can infect up to 10-15 other people through close contact over the course of a year. Due to the highly infectious nature of pulmonary TB, it is important to diagnose and treat the disease very early. Despite the availability of highly effective treatment for decades, TB remains a major global health problem mainly because of poor case detection. The most common method for diagnosing pulmonary TB worldwide is sputum smear microscopy. However sensitivity of direct smear microscopy is low and estimates range from 30% to 70%. It is even lower in case of HIV-infected patients. Culture is more sensitive than microscopy and is considered the current gold standard. Culture requires specialized and controlled laboratory facility and highly skilled manpower and takes 3 to 6 weeks to provide the result. Molecular techniques such as polymerase chain reaction (PCR) or Real Time PCR are much more sensitive than microscopy and culture. However PCR or Real Time PCR tests have so far been restricted to centralized reference laboratories as they require skilled manpower and elaborate infrastructure. Also the turnaround time for results could take a few days.

The **TruelabTM Real Time micro PCR System** enables decentralization and near patient diagnosis of MTB by making real time PCR technology rapid, simple, robust and user friendly and offering "sample to result" capability even at resource limited settings. This is achieved through a combination of lightweight, portable, mains / battery operated **TruelabTM Uno / Uno Dx** Real Time micro PCR Analyzer and **TrueprepTM MAG/AUTO** Sample Prep Device and room temperature stable **TruenatTM micro PCR chip** and **TrueprepTM MAG/AUTO** Sample Prep kits so that even the peripheral laboratories with minimal infrastructure and minimally trained technicians can easily perform these tests routinely in their facilities and report PCR results in less than an hour. Moreover, with these devices PCR testing can also be initiated in the field level, on site.

TruenatTM MTB is a disposable, room temperature stable, chip-based Real Time PCR test with dried down PCR reagents for performing Real Time PCR test for detection of *Mycobacterium tuberculosis* and runs on the **TruelabTM Real Time micro PCR Analyzer**. It requires only six (6) µL of purified DNA to be added to the reaction well for the analysis. The intelligent chip also carries test and batch related information including standard values for quantitation. The **TruenatTM MTB** chip-based Real Time PCR test also stores information of used test to prevent any accidental re-use of the test.



NOTE : **TruelabTM / TruelabTM Uno / TruelabTM Uno Dx / TrueprepTM AUTO/ TrueprepTM MAG / TruepetTM / TruenatTM are all registered trademarks of Molbio Diagnostics (P) Limited.**

The **TruelabTM Real Time micro PCR Analyzer** is protected by the following patents and patents pending: IN 2313/CHE/2007, WO 2009/047804 and corresponding claims of any foreign counterpart(s) thereof. The **TruenatTM micro PCR chip** is protected by the following patents and patents pending: IN 2312/CHE/2007, WO 2009/047805 and corresponding claims of any foreign counterpart(s) thereof. The **TruenatTM MTB chip-based Real Time PCR test** is protected by the following patents and patents pending: IN 796/CHE/2012 and corresponding claims of any foreign counterpart(s) thereof.

3. PRINCIPLE OF THE TEST

TruenatTM MTB works on the principle of Real Time Polymerase Chain Reaction. The DNA from the patient sputum sample is first extracted using **TrueprepTM MAG** Sample Prep Device and **TrueprepTM MAG** sputum Sample Prep kit or using **TrueprepTM AUTO** Universal Cartridge Based Sample Prep Device and **TrueprepTM AUTO** Universal Cartridge Based Sample Prep Kit. The **TruenatTM MTB** chip-based Real Time PCR test is inserted into the **TruelabTM Analyzer** where thermal cycling takes place. Six (6) µL of the purified DNA is dispensed into the reaction well of the **TruenatTM MTB** chip-based Real Time PCR test and the test is started. A positive amplification causes the dual labeled fluorescent probe in the **TruenatTM MTB** chip-based Real Time PCR test to release the fluorophore in an exponential manner which is then captured by the built-in optoelectronic sensor and displayed as amplification curve on the analyzer screen, on a real time basis during the test run. The Cycle threshold (Ct) is defined as the number of amplification cycles required for the fluorescent signal to cross the threshold (i.e. exceed the background signal). Ct levels are inversely proportional to the amount of target nucleic acid in the sample. (i.e. the lower the Ct level the greater is the amount of target nucleic acid in the sample). In the case of negative samples, amplification does not occur and a horizontal amplification curve is displayed on the screen during the test run. At the end of the test run, a MTB "DETECTED" or "NOT DETECTED" result is displayed and in positive cases, quantitative values is also displayed on the screen. Based on the Ct of the internal positive control (IPC), the validity of the test run is also displayed. The IPC is a full process control that undergoes all the processes the specimen undergoes - from extraction to amplification thereby validating the test run from sample to result. Absence of or shift of IPC Ct beyond a pre-set range in case of negative samples invalidates the test run. While IPC will co-amplify in most positive cases also, in some specimen having a high target load, the IPC may not amplify, however the test run is still considered valid. The results can be printed via Bluetooth using the **TruelabTM micro PCR printer** or transferred to the lab computer/ or any remote computer via Wifi network or 3G/GPRS network. Upto 5000 results in **TruelabTM Uno** to 20,000 results in **TruelabTM Uno Dx** can be stored on the analyzer for future recall and reference.

4. TARGET SELECTION

The target sequence for this kit is part of the ribonucleoside-diphosphate reductase gene, the product of which provides the precursor for DNA synthesis. The region selected is specific to the MTB complex.

5. CONTENTS OF THE TruenatTM MTB KIT

A. Individually sealed pouches, each containing

1. **TruenatTM MTB** micro PCR chip.
2. DNase & RNase free pipette tip.
3. Desiccant pouch.

B. Package Insert.

REF	601030005	601030020
	5T	20T

6. CONTENTS OF TrueprepTM AUTO MTB Sample pre-treatment pack (only for TrueprepTM AUTO users)

- A. Liquefaction buffer
- B. Lysis buffer

C. Disposable transfer pipette (graduated)

REF	60204AS05	60204AS20
	5T	20T

7. STORAGE AND STABILITY

TruenatTM MTB chip is stable for one year from the date of manufacture if stored between 2-30°C. It is also stable for upto one (1) month at temperatures up to 40°C. Avoid exposure to light or elevated temperatures (above recommended levels).

TrueprepTM AUTO MTB Sample pretreatment pack is stable for two years from the date of manufacture if stored between 2-30°C. It is also stable for 4 weeks at temperatures up to 45°C. Do not freeze.

8. MATERIALS REQUIRED BUT NOT PROVIDED WITH THE KIT

TruelabTM Real Time micro PCR Workstation (REF 603010001/613010001/623010001) consisting of

1. **TrueprepTM MAG/ TrueprepTM AUTO** Sample Prep Device (REF 603040001/ 603041001).
2. **TruelabTM Uno/ TruelabTM Uno Dx** Real Time micro PCR Analyzer (REF 603020001/ 603021001).
3. **TruelabTM micro PCR Printer** (REF 603050001).
4. **TruepetTM Precision Micropipettes**.

Also required additionally are: **TrueprepTM MAG Sputum Sample Prep Kit** (REF 602020050)/ **TrueprepTM AUTO Universal Cartridge Based Sample Prep Kit** (REF 60203AR25), **TruenatTM Universal Control Kit** (REF 601100008), DNase and RNase-free pipette tips with filter barrier, which may also be procured from **Molbio**, Powder free disposable gloves, waste disposal container with lid.

9. SPECIMEN PREPARATION FOR EXTRACTION WITH TRUEPREPTM MAG

TruenatTM MTB requires purified nucleic acids from pulmonary (sputum/non-sputum) and EPTB specimen that are extracted using the **TrueprepTM MAG** Sample Prep Device and **TrueprepTM MAG** Sputum Sample Prep Kit (Refer to the User Manual of **TrueprepTM MAG** Sample Prep Device and the package insert of **TrueprepTM MAG** Sputum Sample Prep Kit for details).

10. SPECIMEN PREPARATION FOR EXTRACTION WITH TRUEPREPTM AUTO

TruenatTM MTB requires purified nucleic acids from pulmonary (sputum/non-sputum) and EPTB specimen that are extracted using the **TrueprepTM AUTO** Universal Cartridge Based Sample Prep Device and **TrueprepTM AUTO** Universal Cartridge Based Sample Prep Kit. Samples must be liquefied and pre-treated using the **TrueprepTM AUTO** MTB sample pre-treatment pack provided, as per protocol below, before proceeding for extraction.

For sputum samples

Check if the specimen is pipettable. If not, add 1 drop of liquefaction buffer to the specimen (If specimen is frozen allow it to reach room temperature first). Allow the reagent to hydrate the sample by swirling gently. Incubate at room temperature for 5 minutes. If sample has not liquefied after 5 minutes, incubate for another 5 minutes until sample is pipettable. This depends on sample viscosity and ambient temperature.

Label a lysis buffer bottle with patient ID and transfer 500 µl of the liquefied sample into the lysis buffer bottle using the graduated disposable transfer pipette provided. Add 2 drops of the liquefaction buffer to the lysis buffer bottle and mix gently. Close the cap tightly and mix well. Wait for 3 minutes.

Check if contents are fully liquefied by shaking the bottle. If not, incubate it further till the contents are liquefied. Depending upon the sample, this may take another 10 to 15 minutes. Do not proceed if the content has not liquefied.

Sample Storage and Transportation:

Sample pre-treatment decontaminates the specimen and makes it ready for extraction. Sample in this form is stable for 3 days at upto 40°C, and 1 week at 30°C.

Nucleic acid extraction: Follow Extraction procedure (section-13) of **TrueprepTM AUTO** Universal Cartridge Based Sample Prep Kit package insert. (Refer to the User Manual of **TrueprepTM AUTO** Universal Cartridge Based Sample Prep Device and the package insert of **TrueprepTM AUTO** Universal Cartridge Based Sample Prep Kit for details).

For non-sputum samples

Refer protocol sheet in the **TrueprepTM AUTO** MTB Sample Pre-Treatment pack.

Nucleic acid extraction: Follow Extraction procedure (section-13) of **TrueprepTM AUTO** Universal Cartridge Based Sample Prep Kit package insert. (Refer to the User Manual of **TrueprepTM AUTO** Universal Cartridge Based Sample Prep Device and the package insert of **TrueprepTM AUTO** Universal Cartridge Based Sample Prep Kit for details).

11. SAFETY PRECAUTIONS

1. For *in vitro* diagnostic use only.
2. Bring all reagents and specimen to room temperature (20 - 30°C) before use.
3. Do not use kit beyond expiry date.
4. Carefully read the User Manuals and package inserts of all the components of the **TruelabTM Real Time micro PCR System** before use.
5. All materials of human origin should be handled as though potentially infectious.
6. Do not pipette any material by mouth.
7. Do not eat, drink, smoke, apply cosmetics or handle contact lenses in the area where testing is done.
8. Use protective clothing and wear disposable gloves when handling samples and while performing sample extraction.

12. PROCEDURAL PRECAUTIONS

1. Check all packages before using the kit. Damage to the packaging does not prevent the contents of the kit from being used. However, if the outer packaging is damaged the user must confirm that individual components of the kit are intact before using them.
2. Do not perform the test in the presence of reactive vapours (e.g., from sodium hypochlorite, acids, alkalis or aldehydes) or dust.
3. While retrieving the **TruenatTM MTB** micro PCR chip and the DNase & RNase free pipette tip from the pouch, ensure that neither bare hands nor gloves that have been used for previous tests run are used.

13. PROCEDURAL LIMITATIONS

1. Optimal performance of this test requires appropriate specimen collection, handling, storage and transport to the test site.
2. Though very rare, mutations within the highly conserved regions of the target genome where the **TruenatTM assay primers and/or probe** bind may result in the under-quantitation of or a failure to detect the presence of the concerned pathogen.
3. The instruments and assay procedures are designed to minimize the risk of contamination by PCR amplification products. However, it is essential to follow good laboratory practices and ensure careful adherence to the procedures specified in this package insert for avoiding nucleic acid contamination from previous amplifications, positive controls, or specimens.
4. A specimen for which the **TruenatTM assay** reports "Not Detected" cannot be concluded to be negative for the concerned pathogen. As with any diagnostic test, results from the **TruenatTM assay** should be interpreted in the context of other clinical and laboratory findings.

14. CLEANING AND DECONTAMINATION

1. Spills of potentially infectious material should be cleaned up immediately with absorbent paper tissue and the contaminated area should be decontaminated with disinfectants such as 0.5% freshly prepared sodium hypochlorite (10 times dilution of 5% sodium hypochlorite (household bleach) before continuing work.

- Sodium hypochlorite should not be used on an acid-containing spill unless the spill-area is wiped dry first. Materials used to clean spills, including gloves, should be disposed off as potentially bio-hazardous waste e.g. in a biohazard waste container.

15. TEST PROCEDURE

(Please also refer the **Truelab™ Uno / Truelab™ Uno Dx** Real Time micro PCR Analyzer user manual)

- Switch on the **Truelab™** Analyzer.
- If using the **Truelab™ Uno** device, also switch on the touch screen. If using the **Truelab™ Uno Dx** proceed to step 3.
- Select user and enter password.
- Select the test profile for "MTB" on the Analyzer screen.
- Enter the patient details as prompted in the **Truelab™** Analyzer screen.
- Press Start Reaction.
- Press the eject button to open the chip tray.
- Open a pouch of **Truenat™ MTB** and retrieve the chip-based Real Time PCR test.
- Label the chip with the patient ID using a marker pen at the space provided on the back side of the chip.
- Place the **Truenat™ MTB** chip on the chip tray without touching the white reaction well. The reaction well should be facing up and away from the Analyzer. Gently press the chip to ensure that it is seated in the chip tray properly.
- Using the filter barrier tip provided in the pouch, pipette six (6) µL of the purified DNA from the Elute Collection Tube into the centre of the white reaction well. Take care not to scratch the internal well surface and not to spill elute on the outside of the well.
- Slide the chip tray containing the **Truenat™ MTB** chip-based Real Time PCR test loaded with the sample, into the **Truelab™** Analyzer.
- Press **Done** on the "Please Load Sample" Alert message.
- Read the result from the screen.
- Take out the **Truenat™ MTB** micro PCR chip at end of the test and dispose it off as per the section on "Disposal and Destruction" (Section 18).
- Turn on **Truelab™** micro PCR printer and select print on the screen for printing out hard copy of the results. Test results are automatically stored and can be retrieved any time later. (Refer to **Truelab™** Analyzer manual).
- Switch off the **Truelab™** Analyzer.

16. RESULTS & INTERPRETATIONS

Two amplification curves are displayed on the **Truelab™** Analyzer screen when optical plot is selected to indicate the progress of the test. Both the target and the internal positive control (IPC)* curves will take a steep, exponential path when the fluorescence crosses the threshold value in case of positive samples. The Ct will depend on the number of bacterial genomes in the sample. The target curve will remain horizontal throughout the test duration and the IPC curve will take an exponential path in case of negative samples. In case the IPC curve remains horizontal in a negative sample, the test is considered as Invalid. At the end of the test run, the results screen will display "DETECTED" for Positive result or "NOT DETECTED" for Negative result. The result screen would also display the Ct value and the colony forming units per milliliter (cfu/ml) for positive specimen. The result screen also displays the validity of the test run as "VALID" or "INVALID". Invalid samples have to be repeated with fresh specimen from the sample preparation stage. *Note: IPC will co-amplify in most positive cases also, in some specimen having a high target load, the IPC may not amplify, however the test run is still considered valid.

17. QUALITY CONTROL PROCEDURES

To ensure that the **Truelab™** Analyzer is working accurately, run positive and negative controls from time to time. The **Truenat™** Universal Control Kit containing Positive Control and Negative Control must be ordered separately. It is advisable to run controls under the following circumstances:

- Whenever a new shipment of test kits is received.
- When opening a new test kit lot.
- If the temperature of the storage area falls outside of 2-30° C.
- By each new user prior to performing testing on clinical specimen.

18. DISPOSAL AND DESTRUCTION

- Submerge the used **Truenat™ MTB** micro PCR chip in freshly prepared 0.5 % sodium hypochlorite solution for 30 minutes before disposal as per the standard medical waste disposal guidelines.
- Disinfect the solutions and/or solid waste containing biological samples before discarding them according to local regulations.
- Samples and reagents of human and animal origin, as well as contaminated materials, disposables, neutralized acids and other waste materials must be discarded according to local regulations after decontamination by immersion in a freshly prepared 0.5% of sodium hypochlorite for 30 minutes (1 volume of 5% sodium hypochlorite for 10 volumes of contaminated fluid or water).
- Do not autoclave materials or solutions containing sodium hypochlorite.
- Chemicals should be handled in accordance with Good Laboratory Practice and disposed off according to the local regulations.

19. SPECIFIC PERFORMANCE CHARACTERISTICS

ASSAY RANGE AND LIMIT OF DETECTION

To determine the efficiency of the PCR master mix and primer/probe set for MTB DNA amplification, standard curves were developed using serial log dilutions of target DNA. Plasmids carrying cloned MTB PCR amplicons were used as template DNA. The plasmids were quantified by UV spectrophotometry before dilutions, to determine the colony forming units per milliliter (cfu/ml). The **Truenat™ MTB** assay is linear over 7 log MTB DNA dilutions and can detect as low as 400 cfu/ml of sputum.

ANALYTICAL INCLUSIVITY USING WHO-TDR Strain Bank (Sensitivity):

All 235 strains in the TDR-TB-Strain Bank were tested in a blinded study to check the sensitivity of **Truenat™ MTB** against a panel of diverse MTB strains.

Specimen used	Provided by	MTB Positive	MTB Negative	Truenat™ MTB result
235 strains of heat inactivated Mycobacterial cell suspension	ITM-Antwerp	233	02	Identified all samples correctly

Truenat™ MTB identified all 233 *M. tuberculosis* strains and obtained a negative result for the 2 non-tuberculous mycobacteria that were included (*M. terrae* and *M. avium*) in the test panel. **Truenat™ MTB** correctly identified MTB DNA in the TDR-TB-Strain-Bank panel.

ANALYTICAL EXCLUSIVITY (Primer Specificity):

Various strains of Non-tuberculous Mycobacteria (NTM) including *Mycobacterium fortuitum*, *Mycobacterium abscessus*, *Mycobacterium kansasii* I, *Mycobacterium avium*, *Mycobacterium goodii*, *Mycobacterium smegmatis* and *Mycobacterium szulgai* were spiked into negative sputum sample and tested using **Truenat™ MTB** chip-based Real Time PCR test. Cultures of following common bacterial isolates: *Klebsiella pneumoniae*, *Streptococcus sp.*, *Acinetobacter sp.*, *Escherichia coli*, *Salmonella typhi*, *Enterococcus sp.*, and *Pseudomonas aeruginosa* were also tested using the **Truenat™ MTB** chip-based Real Time PCR test. The NTM panel and the bacterial isolates were not detected by the **Truenat™ MTB** chip-based Real Time PCR test indicating its specificity for the MTB complex.

CLINICAL VALIDATION

A pilot study was conducted at P. D. Hinduja National Hospital and Medical Research Centre (Nikam, Chaitali,

et al. " Rapid diagnosis of *Mycobacterium tuberculosis* with **Truenat™ MTB**: a near-care approach." PLOS One 8.1(2013): e51121). 226 sputum specimens from suspected TB patients were analyzed using smear microscopy, culture, in-house nested PCR and **Truenat™ MTB**. Pelleted sputum specimens were re-suspended in lysis buffer from the **Trueprep™ MAG Sputum** kit and processed using the **Trueprep™ MAG Sample Prep Device** followed by PCR on **Truenat™ MTB** chip-based Real Time PCR test. Results were compared with a Composite Reference Standard (CRS) comprising microbiological tests, clinical and radiological findings and patient history. The results are tabulated below:

	Smear		Culture		In-house Nested PCR		Truenat™ MTB	
	+ve	-ve	+ve	-ve	+ve	-ve	+ve	-ve
CRS +ve	120	71	141	50	173	18	174	17
CRS -ve	00	35	00	35	03	32	00	35
Sensitivity %	62.83		73.82		90.58		91.10	
Specificity %	100		100		91.43		100	
PPV%	100		100		98.30		100	
NPV%	33.02		41.18		64.00		67.31	

CRS - Composite Reference Standard, PPV - Positive Predictive Value, NPV - Negative Predictive Value. The results show that **Truenat™ MTB** was the most sensitive (91.10%) and specific (100%) test compared with the Composite Reference Standard. **Truenat™ MTB** also showed high sensitivities of 99.12% among smear positive and culture positive specimen and 75.86% among smear negative and culture positive specimen.

Another study evaluating the **Truenat™ MTB** test was performed using a characterized 100 sample panel from suspected TB patients referred to a hospital in South East Asia. The study involved processing of 500µl of each sputum specimen using the **Trueprep™ MAG Sputum Sample Prep Kit** on the **Trueprep™ MAG Sample Prep Device**. The purified nucleic acids were tested using **Truenat™ MTB** chip-based Real Time PCR test and MTB specific primers and probe on a commercial real-time PCR machine.

Sample Type	Commercial real-time PCR machine result	Truenat™ MTB
S+C+	40/40 (100%)	40/40 (100%)
S-C+	30/40 (75%)	30/40 (75%)
S-C-	0/20 (nil detected)	0/20 (nil detected)

The **Truenat™ MTB** chip-based Real Time PCR test was able to detect 100% of the S+C+ samples (40/40), 75% of S-C+ samples (30/40) and gave a negative result for 100% of the S-C- samples (20/20).

REPRODUCIBILITY AND PRECISION

Assay reproducibility and precision was determined across different users, analyzers and reagent lots. Serial log dilutions of MTB genomic DNA were used for this purpose. The same sample panel was provided to 3 users, who ran it on 3 different **Truelab™** analyzers with 3 different lots of **Truenat™ MTB** chip-based Real Time PCR test. High inter-user, inter-analyzer and inter-lot reproducibility was observed (avg. standard deviation = 0.3 Ct) and no significant difference was observed in the Ct values obtained from different users, analyzers or lots.

INTERFERING SUBSTANCES

In order to study if there is any inhibition in the amplification due to potential interfering substances, sputum samples containing a low load of MTB cells (as characterized by smear microscopy and quantitative PCR) were spiked with blood and pus upto 5% of the final sample concentration. The DNA was extracted using **Trueprep™ MAG Sputum Sample Prep Kit** and Real Time PCR analysis was performed on the **Truenat™ MTB**. The tests were carried out in duplicates.

Effect of interferences on detection of MTB using **Truenat™ MTB**.

Sample ID	Ct of Sputum spiked with blood Average Ct values	Ct of Sputum spiked with pus Average Ct values	Ct of Control sample without interfering substances
Sample 1	33.95	35.60	34.6
Sample 2	32.55	33.70	33.8
Sample 3	33.70	35.75	35.8
Sample 4	35.35	34.10	35.1

No specific inhibition was observed due to blood or pus, in the Real Time PCR result using **Truenat™ MTB**.

20. REFERENCES

- WHO Fact sheet March 2012. <http://www.who.int/mediacentre/factsheets/fs104/en/>.
- Todar's Online Textbook of Bacteriology - Kenneth Todar, Ph.D.
- WHO report 2011 Global Tuberculosis Control.
- Karen R Steingart et al. (2006) Sputum processing methods to improve the sensitivity of smear microscopy for tuberculosis: a systematic review. The Lancet Infectious Diseases - Volume 6, Issue 10, pp 664 - 674.
- P. Fania et al. (2002) Improving Sensitivity of Direct Microscopy for Detection of Acid-Fast Bacilli in Sputum: Use of Chitin in Mucus Digestion, J. Clin Microbiol.; 40(2): 508-511.
- E. Ogbaini-Emovon (2009) Current Trends In The Laboratory Diagnosis of Tuberculosis, Benin Journal of Post Graduate Medicine, Vol.11 Supplemental, pp. 79 - 90.
- Dye C, Watt C.J., Bleed D.M. et al. (2005). Evolution of tuberculosis control and prospects for reducing tuberculosis incidence, prevalence, and deaths globally. JAMA 2005; 293: 2767-2775.
- N. Kennedy et al. (1994) Polymerase Chain Reaction for Assessing Treatment Response in Patients with Pulmonary Tuberculosis. The Journal of Infectious Diseases Vol. 170, No. 3, pp. 713-716.
- Nikam et al. (in press). PLOS ONE.

SYMBOL KEYS

Consult instructions for use	In vitro Diagnostic Test. Not for medicinal use.	Temperature Limitation	Catalogue Number	For single use only
Manufacturer	Date of Manufacture	Date of Expiry	Batch Number / Lot Number	Contains sufficient for $n-1$ tests
EC REP		Authorised Representative in the European Community		

Molbio Diagnostics Pvt. Ltd.

Registered Office:
H.No. 13, Sagar Society, Dona Paula
Panaji, North Goa, Goa - 403 722, India.
www.molbiodiagnostics.com

Manufacturing Unit:
Plot Nos. M 46-47, Phase III B,
Verna Industrial Estate, Verna,
Goa - 403 722, India.
Email: sales@molbiodiagnostics.com

EC REP

Qarad b.v.b.a. Ciplastraat 3, B-2440 Geel, Belgium