



## Technical Specifications

# Trueprep<sup>TM</sup> AUTO

### Universal Cartridge based Sample Prep Device

Trueprep <sup>TM</sup> AUTO Device	
Principle	Proprietary matrix-based extraction
Operation	Fully Automatic
Display Screen	2 line alphanumeric LCD
Power	Rechargeable Lithium Ion Battery Pack 7.4 V, 11.6 Ah External AC/DC adaptor: Input 1.5A, 100/240V, 50/60 Hz; Output 10V, 4.5A
Weight	2.9 kgs
Size	215 x 235 x 115 mm
Software	Proprietary firmware
Operating environment	Temperature 15 - 45°C, RH: 10 -90%

**The True Game Changer**